## **ANIMAL HYPOXIC GLVOE BOX**



## 01 ACCURACY

- Precise O<sub>2</sub> control
  - √ Range of O<sub>2</sub> Control: 0.5% 20.9%, in 0.1% increments
- Precise temperature control (Optional)
  - √ Range of Temperature Control: 5°C above ambient up to 45°C, in 0.1°C increments
- Advanced Airway Design
  - √ Ensures a high level of homogeneity in the workstation's gas distribution



## **02 CONVENIENCE**

- Easy Access Operation System
  - √ The sleeves are made of a soft, skin-friendly material and can access the workstation easily and quickly without vacuum

- Optional Direct-access operation system
  - √ The sleeves are made of a soft, skin-friendly material and an efficient vacuum pump direct access to the workstation only 8-12 seconds.
  - √ Can be operated directly with lab gloves.
  - √ Glove port caps can be "parked" in special holders and do not take up valuable working space.





**Direct-access Operation System** 


 Interlock can be quickly purged by simply pressing the yellow button



 An optional built-in air pump and filter are used to provide the oxygen supply, avoiding the need for a compressed air cylinder

## **03 CUSTOMISED PRODUCT**





Email: sales@maworde.com

www.maworde.com

